

PACKING INSTRUCTIONS
(PRE-ASSEMBLED KIT)
“INFEKTA PAK 1”

UN 4G/CLASS 6.2/05
USA/+AX4770

IMPORTANT WARNING: This packing instruction details correct assembly of this packaging design. Deviation may result in non-compliance and voiding the certification. Infekta is not responsible for any packaging using substitute components or assembled incorrectly.

SHIPPERS RESPONSIBILITY: includes ensuring that inner containers and cushion/absorbent materials are compatible with the contents, and for AIR SHIPMENTS, that the inner containers for liquids are capable of meeting a minimal internal pressure differential of 95 KPA without leakage.

PACKAGING COMPONENTS:

Printed Outer Box (4-1/4” x 4-1/4” x 7-1/4”)
Corrugated “B” perimeter Insert
A 95 KPA, 1-pint (2-1/2” x 6-1/4”) HDPE Screw-Cap Bottle w/Lined Cap
5-Hole Foam Insert (OPTIONAL)
A White fiberboard Screw-Cap mailing tube, seamed bottom & metal lid (3-1/8” x 6-3/8”)

ASSEMBLY INSTRUCTIONS:

1. Remove Mailing Tube/Plastic Bottle from the kit.
2. Remove the plastic bottle from inside the mailing tube
3. Remove cap from plastic bottle.
4. Place primary vials in cushioning material into the plastic bottle.
Insert absorbent with the primaries
(Optional) Insert 5-hole foam insert into the plastic bottle.
Place primaries in the cutouts of the foam insert-add absorbent.
5. Tightly secure the lined cap onto the plastic bottle (70 lbs torque recommended)
Adding a strip of tape around the cap is suggested to insure security of the specimen.
6. Place the plastic bottle inside the mailing tube & secure the metal lid on the tube
7. Insert the mailing tube w/bottle into the preprinted Infekta Pak 1 box.
Secure the top & bottom flaps of the box with a sealing tape.
8. A SHIPPERS DECLARATION FORM IS REQUIRED WITH THIS PACKAGE!
Attached the form on the blank side of the box with an adhesive envelope.