

PACKING INSTRUCTIONS
(PRE-ASSEMBLED KIT)
“INFEKTA PAK 3”

UN 4G/CLASS 6.2/06
USA/+AX4932

IMPORTANT WARNING: This packing instruction details correct assembly of this packaging design. Deviation may result in non-compliance and voiding the certification. Infekta is not responsible for any packaging using substitute components or assembled incorrectly.

SHIPPERS RESPONSIBILITY: includes ensuring that inner containers and cushion/absorbent materials are compatible with the contents, and for AIR SHIPMENTS, that the inner containers for liquids are capable of meeting a minimal internal pressure differential of 95 KPA without leakage.

PACKAGING COMPONENTS:

- Printed Outer Box (4-1/4” x 4-1/4” x 11-1/4”)
- Corrugated “B” flute perimeter Insert
- A 95 KPA, 7 mil Absorbent Lined Specimen Bag (7-1/2” x 11”), printed w/Biohazard Symbol
- A Brown fiberboard Screw-Cap mailing tube, seamed bottom & metal lid (3-1/8” x 11-1/8”)

ASSEMBLY INSTRUCTIONS:

1. Remove Mailing Tube/Specimen Bag from the kit.
2. Remove the specimen bag from inside the mailing tube
3. Place primary vials in cushioning material into the specimen bag.
With the bag flap reverse folded, roll your fingers upward to release the adhesive liner.
4. Lay the bag on a flat surface, fold the bag flap over onto the body of the bag to seal it.
5. Be sure to firmly press down on the entire flap area to adequately seal the adhesive flap to the body of the bag.
6. Form a semi-circle, roll the bag, to accommodate insertion into the mailing tube. Place it into the mailing tube and secure the metal lid of the tube
7. Insert the mailing tube into the preprinted Infekta Pak 3 box.
Secure the top & bottom flaps of the box with a sealing tape.
8. A SHIPPERS DECLARATION FORM IS REQUIRED WITH THIS PACKAGE!
Attached the form on the blank side of the box with an adhesive envelope.