

PACKING INSTRUCTIONS
(PRE-ASSEMBLED KIT)
“INFEKTA PAK 1 (4 PACK)”

UN 4G/CLASS 6.2/05
USA/+AX4771

IMPORTANT WARNING: This packing instruction details correct assembly of this packaging design. Deviation may result in non-compliance and voiding the certification. Infekta is not responsible for any packaging using substitute components or assembled incorrectly.

SHIPPERS RESPONSIBILITY: includes ensuring that inner containers and cushion/absorbent materials are compatible with the contents, and for AIR SHIPMENTS, that the inner containers for liquids are capable of meeting a minimal internal pressure differential of 95 KPA without leakage.

PACKAGING COMPONENTS:

- Printed Outer Box (8-1/2” x 8-1/2” x 7-1/4”)
- Corrugated 4 Chamber Divider/Insert
- (4) 95 KPA, 1-pint (2-1/2” x 6-1/4”) HDPE Screw-Cap Bottles w/Lined Caps
- 5-Hole Foam Insert (OPTIONAL)
- (4) White fiberboard Screw-Cap mailing tubes, seamed bottom & metal lids (3-1/8” x 6-3/8”)

ASSEMBLY INSTRUCTIONS:

1. Remove Mailing Tubes/Plastic Bottles from the kit.
2. Remove the plastic bottles from inside the mailing tubes
3. Remove caps from plastic bottles.
4. Place primary vials in cushioning material into the plastic bottles.
Insert absorbent with the primaries
(Optional) Insert 5-hole foam inserts into the plastic bottles.
Place primaries in the cutouts of the foam inserts-add absorbent.
5. Tightly secure the lined caps onto the plastic bottles (70 lbs torque recommended).
Adding a strip of tape around the caps is suggested to insure security of the specimen.
6. Place the plastic bottles inside the mailing tubes & secure the metal lids on the tubes.
7. Insert the mailing tube w/bottle into each chamber of the preprinted Infekta Pak 1(4 Pack) box.
Secure the top flap of the box with a sealing tape.
8. A SHIPPERS DECLARATION FORM IS REQUIRED WITH THIS PACKAGE!
Attached the form on the blank side of the box with an adhesive envelope.