

## PACKING INSTRUCTIONS BIOLOGICAL CATEGORY B SPECIMENS

**IMPORTANT WARNING:** This packing instruction details correct assembly of this packaging design. Deviation may result in non-compliance and voiding of the certification. ICC The Compliance Center Inc. is not responsible for any packagings using substitute components or assembled incorrectly. If further information is required, please contact ICC The Compliance Center Inc. at 1-800-767-7231

**SHIPPERS RESPONSIBILITY:** includes insuring that all inner containers, closure, liners, cushioning, and (absorbent when applicable) will not react dangerously with the contents, and for AIR SHIPMENT, that inner containers for liquids are capable of meeting a minimum internal pressure differential of 95 kpa without leakage. Some volatile liquids may require a higher test pressure. See ICAO Technical Instruction, Part 3, Chapter 1, Sec 1.1.6.1 or IATA DGR 5.0.2.9

### PACKAGING COMPONENTS:

1. White Fiberboard Box-BX-106A 11 1/4"x 6 5/8"x 3 1/2"
2. Primary Bag (9"x 12") with Biohazard Marking
3. Secondary Bag-95 kPa self sealing bag (6"x 8 1/2")
4. Absorbent Strip 4 1/2"x 2 1/4"



### Packaging Assembly Instructions

1. Place specimens in Primary Bag along with absorbent strip and seal bag.
2. Insert the Primary Bag into the 95 kPa envelope and follow closing instructions.
3. Write any contents and any special instructions on 95 kPa envelope.
4. Place entire envelope into White Box and close box.
5. Follow box closing instructions illustrated below.(1-7)

1



2



3




4



5



6



7

